

STERILE CONTAINER SYSTEMS



A1 Model

STERILE CONTAINER SYSTEMS

MEDIKA
Medizintechnik GmbH



Seven
Different
Models

A1 MODEL



SUITABILITY OF OUR STERILE CONTAINERS FOR STEAM STERILIZATION PROCEDURES ¹⁾
as per DIN 58953 part 9

Sterilization container	Packing		Steam sterilization procedure acc. to DIN 58946, part 1				
	Description	Symbolic representation and packing system	Gravitation procedure	Pre-vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
with filters (StB-F)	Textiles in sterilization container with filter in lid and bottom acc. to DIN 58952 part 1		+	+	+	+	+
			○	○	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated acc.to DIN 58952 part 1		-	-	+	+	+
			-	-	+	+	○
	Textiles in sterilization container with filter in bottom lid, not perforated acc.to DIN 58952 part 1		○	+	+	+	+
			-	○	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid and bottom acc. to DIN 58952 part 1		+ ³⁾	+	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid and bottom not perforated acc. to DIN 58952 part 1		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in bottom acc. to DIN 58952 part 1		○	+	+	+	+
with Valves (StB-V)	Textiles in sterilization container with valves in lid and bottom acc. to DIN 58952 part 1		-	○	+	+	○
			-	-	+	+	○
	Textiles in sterilization container with valves in lid no valves on the bottom		-	-	+	+	-
			-	-	+	+	-
	Instruments in wire basket in sterilization container with valve in lid and bottom		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket in sterilization container with valve in lid and not perforated		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾

Explanation of symbols:

+ applicable

- not applicable

○ limited application

1) In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.

2) Only sterilization containers guaranteeing evaporation or discharge of condensate are suitable above.

3) Gravitation procedures having no drying device do not guarantee a sufficient after-drying effect.




STERILE CONTAINER SYSTEMS









A1 Model

» 1/1 Full Size

1/1 FULL SIZE



h=100 mm	 K6-010E	 K6-011E	 K6-012E	 K6-013E	 K6-014E
h=135 mm	 K6-016E	 K6-017E	 K6-018E	 K6-019E	 K6-020E
h=150 mm	 K6-022E	 K6-023E	 K6-024E	 K6-025E	 K6-026E
h=200 mm	 K6-028E	 K6-029E	 K6-030E	 K6-031E	 K6-032E


BOX (non perforated)


outside dimension	inside dimension	weight
600 x 272 h=100 mm	543 x 258 h=94 mm	2.00 kg
600 x 272 h=116 mm	543 x 258 h=112 mm	2.10 kg
600 x 272 h=138 mm	543 x 258 h=130 mm	2.32 kg
600 x 272 h=190 mm	543 x 258 h=180 mm	2.66 kg

PERFORATED LID


Dimension: 602 x 285 mm
Weight : 1.35 kg

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
	K3-740	535 x 240 h=35 mm
	K3-741	535 x 240 h=55 mm
	K3-742	535 x 240 h=75 mm
	K3-743	535 x 240 h=105 mm

STAND



code	height
K3-767	20 mm



STAND



code	height
K3-768	10 mm



STAND



code	height
K3-766	5 mm




BASKET MADE FROM STAINLESS STEEL WIRE MESH




code	outside dimension
K2-460	540 x 250 h=30 mm
K2-461	540 x 250 h=50 mm
K2-462	540 x 250 h=70 mm
K2-463	540 x 250 h=100 mm

SILICONE MAT



K2-526
520 x 230 mm
K2-528
470 x 230 mm

SILICONE MESH MAT



K2-531
480 x 250 mm

SILICONE HANDLE



K2-710	Red
K2-711	Yellow
K2-712	Green
K2-713	Blue
K2-714	Grey
K2-715	Black

PAPER FILTER



K2-740
Single use 100 pcs
with indicator. Aygun mark.
Diameter: 185 mm

ALUMINIUM FILTER RETAINER




K2-722
Diameter: 189 mm

FILTER DISPENSER



K2-730

IDENTIFICATION LABELS (aluminium)



K2-700	Red
K2-701	Blue
K2-702	Green
K2-703	Yellow
K2-704	Grey
K2-705	Black

LABELS (paper)



K2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



K2-672
100 pcs. With indicator.

For full and
Half size containers

SECURITY SEAL (blue)



K2-663
100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)



K2-664
100 pcs. made of plastic

SECURITY SEAL (blue)



K2-661
100 pcs. with number.
made of plastic

SECURITY SEAL (red)



K2-662
100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)



K2-660
100 pcs. made of plastic

SECURITY SEAL (red)



K2-665
100 pcs. with number.
made of plastic

A1 MODEL

1/1 full size






STERILE CONTAINER SYSTEMS





A1 Model

» 3/4 Quarter Size


3/4 QUARTER SIZE




h=100 mm	 K6-240E	 K6-241E	 K6-242E	 K6-243E	 K6-244E
h=135 mm	 K6-246E	 K6-247E	 K6-248E	 K6-249E	 K6-250E
h=150 mm	 K6-252E	 K6-253E	 K6-254E	 K6-255E	 K6-256E


BOX (non perforated)



outside dimension



outside dimension



inside dimension



outside dimension	inside dimension	weight
485 x 272 h=100 mm	422 x 258 h=94 mm	1.64 kg
485 x 272 h=116 mm	422 x 258 h=112 mm	1.78 kg
485 x 272 h=138 mm	422 x 258 h=133 mm	1.90 kg

PERFORATED LID


Dimension: 477 x 285 mm
Weight : 1.05 kg

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
	K3-748	400 x 240 h=35 mm
	K3-749	400 x 240 h=55 mm
	K3-750	400 x 240 h=75 mm
	K3-751	400 x 240 h=105 mm


STAND



code	height
K3-767	20 mm



STAND



code	height
K3-768	10 mm



STAND



code	height
K3-766	5 mm



BASKET MADE FROM STAINLESS STEEL WIRE MESH




code	outside dimension
K2-1500	405 x 250 x h=30 mm
K2-1501	405 x 250 x h=50 mm
K2-1502	405 x 250 x h=70 mm
K2-1503	405 x 250 x h=100 mm

SILICONE MAT



K2-524
 380 x 230 mm
K2-525
 440 x 230 mm

SILICONE MESH MAT



K2-531
 480 x 250 mm

SILICONE HANDLE



K2-710	
K2-711	
K2-712	
K2-713	
K2-714	
K2-715	

PAPER FILTER



K2-740
 Single use 100 pcs
 with indicator. Aygun mark.
 • Diameter: 185 mm

ALUMINIUM FILTER RETAINER




K2-722
 Diameter: 189 mm

FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)



K2-700	
K2-701	
K2-702	
K2-703	
K2-704	
K2-705	

LABELS (paper)



K2-673
 100 pcs. With indicator.

For small size and
 special containers

LABELS (paper)



K2-672
 100 pcs. With indicator.

For full and
 Half size containers

SECURITY SEAL (blue)



K2-663
 100 pcs. with indicator.
 made of plastic

SECURITY SEAL (blue)



K2-664
 100 pcs. made of plastic

SECURITY SEAL (blue)



K2-661
 100 pcs. with number.
 made of plastic

SECURITY SEAL (red)



K2-662
 100 pcs. with indicator.
 made of plastic

SECURITY SEAL (red)



K2-660
 100 pcs. made of plastic

SECURITY SEAL (red)



K2-665
 100 pcs. with number.
 made of plastic

A1 MODEL

3/4 Quarter size






STERILE CONTAINER SYSTEMS





A1 Model

» *1/2 Half Size*


1/2 HALF SIZE




h=100 mm	 K6-500E	 K6-501E	 K6-502E	 K6-503E	 K6-504E
h=135 mm	 K6-506E	 K6-507E	 K6-508E	 K6-509E	 K6-510E
h=150 mm	 K6-512E	 K6-513E	 K6-514E	 K6-515E	 K6-516E
h=200 mm	 K6-518E	 K6-519E	 K6-520E	 K6-521E	 K6-522E


BOX (non perforated)



outside dimension



outside dimension



inside dimension



outside dimension	inside dimension	weight
305 x 272 h=100 mm	253 x 258 h=94 mm	1.22 kg
305 x 272 h=116 mm	253 x 258 h=112 mm	1.38 kg
305 x 272 h=138 mm	253 x 258 h=133 mm	1.49 kg
305 x 272 h=190 mm	253 x 258 h=184 mm	1.67 kg

PERFORATED LID


Dimension: 308 x 285 mm
Weight : 0.75 kg

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
	K3-752	244 x 235 h=35 mm
	K3-753	244 x 235 h=55 mm
	K3-754	244 x 235 h=75 mm
	K3-755	244 x 235 h=105 mm

STAND



code	height
K3-767	20 mm



STAND



code	height
K3-768	10 mm



STAND



code	height
K3-766	5 mm



BASKET MADE FROM STAINLESS STEEL WIRE MESH




code	outside dimension
K2-490	255 x 250 x h=30 mm
K2-491	255 x 250 x h=50 mm
K2-492	255 x 250 x h=70 mm
K2-493	255 x 250 x h=100 mm

SILICONE MAT



K2-521
250 x 235 mm

SILICONE MESH MAT



K2-530
250 x 250 mm

SILICONE HANDLE



K2-710	Red
K2-711	Yellow
K2-712	Green
K2-713	Blue
K2-714	Grey
K2-715	Black

PAPER FILTER



K2-740
Single use 100 pcs
with indicator. Aygun mark.
• Diameter: 185 mm

ALUMINIUM FILTER RETAINER




K2-722
Diameter: 189 mm

FILTER DISPENSER



K2-730

IDENTIFICATION LABELS (aluminium)



K2-700	Red
K2-701	Blue
K2-702	Green
K2-703	Yellow
K2-704	Grey
K2-705	Black

LABELS (paper)



K2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



K2-672
100 pcs. With indicator.

For full and
Half size containers

SECURITY SEAL (blue)



K2-663
100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)



K2-664
100 pcs. made of plastic

SECURITY SEAL (blue)



K2-661
100 pcs. with number.
made of plastic

SECURITY SEAL (red)



K2-662
100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)



K2-660
100 pcs. made of plastic

SECURITY SEAL (red)



K2-665
100 pcs. with number.
made of plastic

A1 MODEL

1/2 half size

