






STERILE CONTAINER SYSTEMS



Dental, Mini & Special Model

LARGE DENTAL MODEL CONTAINER




Specifications Sizes		Lid (Perforated)	
			
		Box (Non Perforated)	
h=100 mm		K6-1600	
		K6-1601	
		K6-1602	
		K6-1603	
		K6-1604	
		K6-1605	
h=135 mm		K6-1606	
		K6-1607	
		K6-1608	
		K6-1609	
		K6-1610	
		K6-1611	
h=150 mm		K6-1612	
		K6-1613	
		K6-1614	
		K6-1615	
		K6-1616	
		K6-1617	


Ps: All container measurements cover the outer size of the lids.

BOX (non perforated)


outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
K6-1665	380 x 244 h=100 mm	330 x 234 h=94 mm	1,35 kg
K6-1666	380 x 244 h=116 mm	330 x 234 h=112 mm	1,48 kg
K6-1667	380 x 244 h=138 mm	330 x 234 h=130 mm	1,55 kg

PERFORATED LIDS


code	color	weight
K6-1680		0,75 kg
K6-1681		
K6-1682		
K6-1683		
K6-1684		
K6-1685		

LID-BARRIER SYSTEM


code	color	weight
K7-1680		0,95 kg
K7-1681		
K7-1682		
K7-1683		
K7-1684		
K7-1685		

 STANDARD
MODEL
1

 BIO-BARRIER
MODEL
33

 Polymer's PPSU
Lids
65

 WINDOW
MODEL
85

 B₃^{Strirrer}
MODEL
103

 for STERRAD
135

 Dental, Mini,
Special
145

 BASKETS
163

 Autolock Mechanism
Emulating Indicator
183

 ACCESSORIES
191

BASKET MADE FROM STAINLESS STEEL WIRE MESH

WIRE BASKET




code	outside dimension
K3-2000	320 x 220 h=50 mm
K3-2001	320 x 220 h=70 mm
K3-2002	320 x 220 h=100 mm

WIRE BASKET WITH LID



code	outside dimension
K3-2010	320 x 220 h=55 mm
K3-2011	320 x 220 h=75 mm
K3-2012	320 x 220 h=105 mm
K3-2013	ONLY LID, for baskets

SILICONE MAT



K2-532
280 x 200 mm

DISPOSAL STAMP



K2-717
40 x 35 x 9 mm
Requested number or letter must be informed during the order.

LABEL HOLDER




K2-716
50 x 35 x 13 mm
Holder for identification labels

PAPER FILTER



K2-740
Single use 100 pcs with indicator. Aygun mark.
• Diameter: 185 mm

TEFLON FILTER



K2-742
Number of sterilization cycle > 2000
• Diameter: 194 mm

FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)



K2-700 (red)
K2-701 (blue)
K2-702 (green)
K2-703 (yellow)
K2-704 (grey)
K2-705 (black)

LABELS (paper)



K2-673
100 pcs. With indicator.

LABELS (paper)



K2-672
100 pcs. With indicator.

SECURITY SEAL (blue)



K2-663
100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)



K2-664
100 pcs. made of plastic

SECURITY SEAL (blue)



K2-661
100 pcs. with number.
made of plastic

SECURITY SEAL (red)



K2-662
100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)



K2-660
100 pcs. made of plastic

SECURITY SEAL (red)




K2-665
100 pcs. with number.
made of plastic

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
Box and Lid	K3-774	320 x 220 h=63 mm
	K3-775	320 x 220 h=83 mm
	K3-776	320 x 220 h=113 mm
Box	K3-771	320 x 220 h=60 mm
	K3-772	320 x 220 h=80 mm
	K3-773	320 x 220 h=110 mm
Lid	K3-770	ONLY LID, for baskets



STAND



code	height
K3-767	20 mm



STAND



code	height
K3-768	10 mm



IDENTIFICATION (aluminium)

LABEL HOLDER



code	color
K2-560	Red
K2-561	Blue
K2-562	Green
K2-563	Yellow
K2-564	Grey
K2-565	Black

LABELS



code	color
K2-570	Red
K2-571	Blue
K2-572	Green
K2-573	Yellow
K2-574	Grey
K2-575	Black

Coloured identification labels
Engraved to your specification



➤ Grey coloured recommended for use in Sterrad.
Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.

Specifications ➔		Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
Sizes						
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)
320 X 190 X 50 mm 	K6-1800 	K6-1820 	K6-1840 	K6-1860 	K6-1865 	
	K6-1801 	K6-1821 	K6-1841 			
	K6-1802 	K6-1822 	K6-1842 			
	K6-1803 	K6-1823 	K6-1843 			
	K6-1804 	K6-1824 	K6-1844 			
	K6-1805 	K6-1825 	K6-1845 			
320 X 190 X 65 mm 	K6-1806 	K6-1826 	K6-1846 	K6-1861 	K6-1866 	
	K6-1807 	K6-1827 	K6-1847 			
	K6-1808 	K6-1828 	K6-1848 			
	K6-1809 	K6-1829 	K6-1849 			
	K6-1810 	K6-1830 	K6-1850 			
	K6-1811 	K6-1831 	K6-1851 			
320 X 190 X 130 mm 	K6-1812 	K6-1832 	K6-1852 	K6-1862 	K6-1867 	
	K6-1813 	K6-1833 	K6-1853 			
	K6-1814 	K6-1834 	K6-1854 			
	K6-1815 	K6-1835 	K6-1855 			
	K6-1816 	K6-1836 	K6-1856 			
	K6-1817 	K6-1837 	K6-1857 			

➔ Ps: All container measurements cover the outer size of the lids.

WIRE BASKET (with lid)



code	outside dimension
K3-611	270 x 175 x h=45 mm

WIRE BASKET




code	outside dimension
K3-1070	270 x 170 x h=45 mm

STORAGE FRAME



code	outside dimension
K3-609	273 x 180 x 13 mm

SILICONE MAT




code	outside dimension
K2-529	260 x 180 mm

DRAPE RETAINER



code	outside dimension
K3-610	277 x 175 mm

PERFORATED BOTTOMS



code	K6-1860
outside dimension	320 x 186 x h ¹ =43 mm
inside dimension	278 x 184 x h ² =28 mm
weight	0,45 kg




code	K6-1861
outside dimension	320 x 186 x h ¹ =58 mm
inside dimension	278 x 184 x h ² =46 mm
weight	0,50 kg



code	K6-1862
outside dimension	320 x 186 x h ¹ =122 mm
inside dimension	278 x 184 x h ² =110 mm
weight	0,70 kg

SOLID BOTTOMS



code	K6-1865
outside dimension	320 x 186 x h ¹ =43 mm
inside dimension	278 x 184 x h ² =38 mm
weight	0,35 kg




code	K6-1866
outside dimension	320 x 186 x h ¹ =59 mm
inside dimension	278 x 184 x h ² =54 mm
weight	0,40 kg




code	K6-1867
outside dimension	320 x 186 x h ¹ =122 mm
inside dimension	278 x 184 x h ² =118 mm
weight	0,65 kg

FILTERS



code **K3-612**
 Dimension: 165 x 235 mm
 Single use **paper filter** 100 pcs.
 With indicator, Aygun mark.



code **K3-615**
 Dimension: 165 x 235 mm
 Single use **paper filter** 100 pcs.
 With indicator, No marking.



code **K2-743**
 Dimension: 150 x 217 mm
Teflon Filter
 Number of sterilization cycle > 2000

PERFORATED LIDS



code	color	weight
K6-1880		0.60 kg
K6-1881		
K6-1882		
K6-1883		
K6-1884		
K6-1885		

REMOVAL LIDS (non perforated)



code	color	weight
K6-1870		0.47 kg
K6-1871		
K6-1872		
K6-1873		
K6-1874		
K6-1875		

BIO-BARRIER LIDS



code	color	weight
K7-1880		0.90 kg
K7-1881		
K7-1882		
K7-1883		
K7-1884		
K7-1885		

STANDARD
MODEL

1

BIO-BARRIER
MODEL

33

Polymers PPSU
Lids

65

WINDOW
MODEL

85

B₃ Sterilizer
MODEL

103

for STERRAD

135

Dental, Mini,
Special

145

BASKETS

163

Autolock Mechanism
Emulating Indicator

183

ACCESSORIES

191

Specifications ➔		Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
Sizes						
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)
236 X 174 X 60 mm		K6-1500 	K6-1520 	K6-1540 	K6-1560 	K6-1565 
		K6-1501 	K6-1521 	K6-1541 		
		K6-1502 	K6-1522 	K6-1542 		
		K6-1503 	K6-1523 	K6-1543 		
		K6-1504 	K6-1524 	K6-1544 		
		K6-1505 	K6-1525 	K6-1545 		
236 X 174 X 100 mm		K6-1506 	K6-1526 	K6-1546 	K6-1561 	K6-1566 
		K6-1507 	K6-1527 	K6-1547 		
		K6-1508 	K6-1528 	K6-1548 		
		K6-1509 	K6-1529 	K6-1549 		
		K6-1510 	K6-1530 	K6-1550 		
		K6-1511 	K6-1531 	K6-1551 		
236 X 174 X 130 mm		K6-1512 	K6-1532 	K6-1552 	K6-1562 	K6-1567 
		K6-1513 	K6-1533 	K6-1553 		
		K6-1514 	K6-1534 	K6-1554 		
		K6-1515 	K6-1535 	K6-1555 		
		K6-1516 	K6-1536 	K6-1556 		
		K6-1517 	K6-1537 	K6-1557 		

➔ Ps: All container measurements cover the outer size of the lids.

WIRE BASKET (with lid)



code	outside dimension
K2-995	170 x 140 x 45 mm

WIRE BASKET



code	outside dimension
K3-1079	170 x 140 x 45 mm

STORAGE FRAME



code	outside dimension
K3-608	170 x 140 mm

SILICONE MAT




code	outside dimension
K2-519	150 x 130 mm

DRAPE RETAINER



code	outside dimension
K3-607	170 x 140 mm

PERFORATED BOTTOMS



code	K6-1560
outside dimension	230 x 161 x h ¹ =50 mm
inside dimension	190 x 158 x h ² =40 mm
weight	




code	K6-1561
outside dimension	230 x 161 x h ¹ =90 mm
inside dimension	190 x 158 x h ² =80 mm
weight	



code	K6-1562
outside dimension	230 x 161 x h ¹ =120 mm
inside dimension	190 x 158 x h ² =110 mm
weight	

SOLID BOTTOMS



code	K6-1565
outside dimension	230 x 161 x h ¹ =50 mm
inside dimension	190 x 158 x h ² =48 mm
weight	



code	K6-1566
outside dimension	230 x 161 x h ¹ =90 mm
inside dimension	190 x 158 x h ² =88 mm
weight	



code	K6-1567
outside dimension	230 x 161 x h ¹ =120 mm
inside dimension	190 x 158 x h ² =118 mm
weight	

FILTERS



code **K2-781**
 Dimension: 130 x 110 mm
 Single use **paper filter** 100 pcs.
 With indicator, Aygun mark.



code **K2-780**
 Dimension: 130 x 110 mm
 Single use **paper filter** 100 pcs.
 With indicator, No marking.



code **K2-782**
 Dimension: 130 x 110 mm
Teflon Filter
 Number of sterilization
 cycle > 2000

PERFORATED LIDS

code	color	weight
K6-1570		
K6-1571		
K6-1572		
K6-1573		
K6-1574		
K6-1575		



REMOVAL LIDS (non perforated)

code	color	weight
K6-1580		
K6-1581		
K6-1582		
K6-1583		
K6-1584		
K6-1585		

STANDARD
MODEL
1BIO-BARRIER
MODEL
33Polymers PPSU
Lids
65WINDOW
MODEL
85B₃ Steril
MODEL
103for STERRAD
135Dental, Mini,
Special
145BASKETS
163Autolock Mechanism
Emulating Indicator
183ACCESSORIES
191

Specifications ➔		Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
Sizes						
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)
320 X 140 X 40 mm		K6-1900 	K6-1920 	K6-1940 	K6-1960 	K6-1965 
		K6-1901 	K6-1921 	K6-1941 		
		K6-1902 	K6-1922 	K6-1942 		
		K6-1903 	K6-1923 	K6-1943 		
		K6-1904 	K6-1924 	K6-1944 		
		K6-1905 	K6-1925 	K6-1945 		
320 X 140 X 70 mm		K6-1906 	K6-1926 	K6-1946 	K6-1961 	K6-1966 
		K6-1907 	K6-1927 	K6-1947 		
		K6-1908 	K6-1928 	K6-1948 		
		K6-1909 	K6-1929 	K6-1949 		
		K6-1910 	K6-1930 	K6-1950 		
		K6-1911 	K6-1931 	K6-1951 		
320 X 140 X 100 mm		K6-1912 	K6-1932 	K6-1952 	K6-1962 	K6-1967 
		K6-1913 	K6-1933 	K6-1953 		
		K6-1914 	K6-1934 	K6-1954 		
		K6-1915 	K6-1935 	K6-1955 		
		K6-1916 	K6-1936 	K6-1956 		
		K6-1917 	K6-1937 	K6-1957 		

➔ Ps: All container measurements cover the outer size of the lids.

WIRE BASKET (with lid)



code	outside dimension
K2-993	265 x 125 x h=50 mm

WIRE BASKET



code	outside dimension
K3-1073	265 x 120 x h=50 mm

WIRE BASE



code	outside dimension
K2-992	270 x 125 mm

SILICONE MAT




code	outside dimension
K2-520	275 x 125 mm

DRAPE RETAINER



code	outside dimension

PERFORATED BOTTOMS



code	K6-1960
outside dimension	315 x 132 x h¹=32 mm
inside dimension	273 x 128 x h²=18 mm
weight	0,30 kg




code	K6-1961
outside dimension	315 x 132 x h¹=63 mm
inside dimension	273 x 128 x h²=48 mm
weight	0,40 kg



code	K6-1962
outside dimension	315 x 132 x h¹=92 mm
inside dimension	273 x 128 x h²=78 mm
weight	0,50 kg

SOLID BOTTOMS



code	K6-1965
outside dimension	315 x 132 x h¹=32 mm
inside dimension	273 x 128 x h²=30 mm
weight	0,20 kg



code	K6-1966
outside dimension	315 x 132 x h¹=63 mm
inside dimension	273 x 128 x h²=58 mm
weight	0,35 kg




code	K6-1967
outside dimension	315 x 132 x h¹=92 mm
inside dimension	273 x 128 x h²=90 mm
weight	0,45 kg

FILTERS



code **K2-760**
 Dimension: 235 x 120 mm
 Single use **paper filter** 100 pcs.
 With indicator, Aygun mark.



code **K2-761**
 Dimension: 235 x 120 mm
 Single use **paper filter** 100 pcs.
 With indicator, No marking.



code **K2-744**
 Dimension: 215 x 100 mm
Teflon Filter
 Number of sterilization cycle > 2000

PERFORATED LIDS



code	color	weight
K6-1970		0.48 kg
K6-1971		
K6-1972		
K6-1973		
K6-1974		
K6-1975		

REMOVAL LIDS (non perforated)



code	color	weight
K6-1980		0.38 kg
K6-1981		
K6-1982		
K6-1983		
K6-1984		
K6-1985		

BIO-BARRIER LIDS



code	color	weight
K7-1970		0.65 kg
K7-1971		
K7-1972		
K7-1973		
K7-1974		
K7-1975		

STANDARD
MODEL

1

BIO-BARRIER
MODEL

33

Polymers PPSU
Lids

65

WINDOW
MODEL

85

B₃
Model

103

for STERRAD

135

Dental, Mini,
Special

145

BASKETS


163

Autolock Mechanism
Emulating Indicator

183

ACCESSORIES

191

Specifications ➔		Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
						
C Sizes		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)
	285 X 280 X 55 mm	 K6-1710  K6-1711  K6-1712  K6-1713  K6-1714  K6-1715 	K6-1730  K6-1731  K6-1732  K6-1733  K6-1734  K6-1735 	K6-1750  K6-1751  K6-1752  K6-1753  K6-1754  K6-1755 	K6-1770 	K6-1775 
	285 X 280 X 85 mm	 K6-1716  K6-1717  K6-1718  K6-1719  K6-1720  K6-1721 	K6-1736  K6-1737  K6-1738  K6-1739  K6-1740  K6-1741 	K6-1756  K6-1757  K6-1758  K6-1759  K6-1760  K6-1761 	K6-1771 	K6-1776 

➔ Ps: All container measurements cover the outer size of the lids.

WIRE BASKET



code	outside dimension
K2-490	255 x 250 x h=30 mm
K2-491	255 x 250 x h=50 mm

SILICONE MAT



code	outside dimension
K2-521	255 x 255 mm

WIRE BASE



code	outside dimension
K2-482	255 x 255 x 30 mm

PAPER FILTER




code	K2-740
Diameter	185 mm
Single use 100 pcs with indicator. Aygun mark.	

DRAPE RETAINER



code	outside dimension
K2-480	255 x 255 mm

PAPER FILTER



code	K2-741
Diameter	185 mm
Single use 100 pcs with indicator. No marking.	

DRAPE RETAINER




code	outside dimension
K2-481	255 x 255 x 30 mm

TEFLON FILTER



code	K2-742
Diameter	194 mm
Number of sterilization cycle > 2000	

PERFORATED BOTTOMS




code	K6-1770
outside dimension	295 x 270 x h ¹ =60 mm
inside dimension	258 x 260 x h ² =40 mm
weight	0,70 kg



code	K6-1771
outside dimension	295 x 270 x h ¹ =72 mm
inside dimension	258 x 260 x h ² =52 mm
weight	0,90 kg

SOLID BOTTOMS



code	K6-1775
outside dimension	295 x 270 x h ¹ =60 mm
inside dimension	258 x 260 x h ² =56 mm
weight	0,60 kg



code	K6-1776
outside dimension	295 x 270 x h ¹ =72 mm
inside dimension	258 x 260 x h ² =66 mm
weight	0,70 kg

PERFORATED LIDS



code	color	weight
K6-1780		0.71 kg
K6-1781		
K6-1782		
K6-1783		
K6-1784		
K6-1785		

REMOVAL LIDS (non perforated)



code	color	weight
K6-1790		0.58 kg
K6-1791		
K6-1792		
K6-1793		
K6-1794		
K6-1795		

BIO-BARRIER LIDS



code	color	weight
K7-1780		0.85 kg
K7-1781		
K7-1782		
K7-1783		
K7-1784		
K7-1785		

SAFETY LIDS



code	color	weight
K6-1418		0.60 kg
K6-1419		
K6-1420		
K6-1421		
K6-1422		
K6-1423		



BIO-BARRIER SAFETY LID



code	color	weight
K7-1418		0.60 kg
K7-1419		
K7-1420		
K7-1421		
K7-1422		
K7-1423		

STANDARD
MODEL
1BIO-BARRIER
MODEL
33Polymers PPSU
Lids
65WINDOW
MODEL
85B₃
Steril
MODEL
103for STERRAD
135Dental, Mini,
Special
145BASKETS
163Autolock Mechanism
Emulating Indicator
183ACCESSORIES
191

ENDO MODEL CONTAINER



color	code	outside dimension
	K6-1400	540 x 159 h=80 mm
	K6-1406	540 x 159 h=100 mm

color	code	outside dimension
	K6-1401	540 x 159 h=80 mm
	K6-1407	540 x 159 h=100 mm


color	code	outside dimension
	K6-1402	540 x 159 h=80 mm
	K6-1408	540 x 159 h=100 mm

color	code	outside dimension
	K6-1403	540 x 159 h=80 mm
	K6-1409	540 x 159 h=100 mm

color	code	outside dimension
	K6-1404	540 x 159 h=80 mm
	K6-1410	540 x 159 h=100 mm

color	code	outside dimension
	K6-1405	540 x 159 h=80 mm
	K6-1411	540 x 159 h=100 mm

PERFORATED BOTTOMS



code	outside dimension
K6-1100	540 x 145 h=67 mm
K6-1101	540 x 145 h=87 mm

PERFORATED LIDS



code	color
K6-1000	
K6-1001	
K6-1002	
K6-1003	
K6-1004	
K6-1005	

FILTERS



code **K2-760**
 Dimension: 235 x 120 mm
 Single use **paper filter** 100 pcs.
 With indicator, Aygun mark.



code **K2-761**
 Dimension: 235 x 120 mm
 Single use **paper filter** 100 pcs.
 With indicator, No marking.



code **K2-744**
 Dimension: 215 x 100 mm
Teflon Filter
 Number of sterilization
 cycle > 2000

BASKET (with stand)



code	outside dimension
K3-1054	500 x 135 x 50 mm



LABELS (paper)



K2-673
 100 pcs. With indicator.

LABELS (paper)



K2-672
 100 pcs. With indicator.

IDENTIFICATION LABELS (aluminium)



K2-700	
K2-701	
K2-702	
K2-703	
K2-704	
K2-705	

STANDARD
MODEL
1BIO-BARRIER
MODEL
33Polymers PPSU
Lids
65WINDOW
MODEL
85B₃ Steril
MODEL
103for STERRAD
135Dental, Mini,
Special
145BASKETS
163Autolock Mechanism
Emulating Indicator
183ACCESSORIES
191

IMPLANT MODEL CONTAINERS


Specifications		Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)
Sizes				
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)
300 X 150 X 75 mm		K5-010 	K5-022 	K5-034 
		K5-011 	K5-023 	K5-035 
		K5-012 	K5-024 	K5-036 
		K5-013 	K5-025 	K5-037 
		K5-014 	K5-026 	K5-038 
		K5-015 	K5-027 	K5-039 
500 X 170 X 75 mm		K5-016 	K5-028 	K5-040 
		K5-017 	K5-029 	K5-041 
		K5-018 	K5-030 	K5-042 
		K5-019 	K5-031 	K5-043 
		K5-020 	K5-032 	K5-044 
		K5-021 	K5-033 	K5-045 

Non filter



Ps: All container measurements cover the outer size of the lids.

MINI IMPLANT MODEL CONTAINERS

Specifications		Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)
Sizes				
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)
160 X 70 X 60 mm		K5-200 	K5-212 	K5-224 
		K5-201 	K5-213 	K5-225 
		K5-202 	K5-214 	K5-226 
		K5-203 	K5-215 	K5-227 
		K5-204 	K5-216 	K5-228 
		K5-205 	K5-217 	K5-229 
200 X 145 X 60 mm		K5-206 	K5-218 	K5-230 
		K5-207 	K5-219 	K5-231 
		K5-208 	K5-220 	K5-232 
		K5-209 	K5-221 	K5-233 
		K5-210 	K5-222 	K5-234 
		K5-211 	K5-223 	K5-235 

Non filter



SCOPY MODEL CONTAINERS

Specifications Sizes	Lid (Perforated)	
	Box (Non Perforated)	Box (Perforated)
260 X 70 X 65 mm	K5-060 K5-061 K5-062 K5-063 K5-064 K5-065	K5-072 K5-073 K5-074 K5-075 K5-076 K5-077
460 X 80 X 65 mm	K5-066 K5-067 K5-068 K5-069 K5-070 K5-071	K5-078 K5-079 K5-080 K5-081 K5-082 K5-083

Ps: All container measurements cover the outer size of the lids.



Non filter

SCOPY MODEL CONTAINERS (with silicon stand, lock, label)

Specifications Sizes	Lid (Perforated)	
	Box (Non Perforated)	Box (Perforated)
270 X 70 X 65 mm	K5-100 K5-101 K5-102 K5-103 K5-104 K5-105	K5-115 K5-116 K5-117 K5-118 K5-119 K5-120
475 X 80 X 65 mm	K5-106 K5-107 K5-108 K5-109 K5-110 K5-111	K5-121 K5-122 K5-123 K5-124 K5-125 K5-126

Ps: All container measurements cover the outer size of the lids.

STANDARD
MODEL
1BIO-BARRIER
MODEL
33Polymers PPSU
Lids
65WINDOW
MODEL
85B₃
Striker
MODEL
103for STERRAD
135Dental, Mini,
Special
145BASKETS
163Autolock Mechanism
Emulating Indicator
183ACCESSORIES
191

BRUSH MODEL CONTAINERS



Container for brush 280 x 285 x 85 mm
With wall mounting fixtures,
Maximum 20 brushes can be stored in container
In order to take the brush, please move the handle under
the container to right or left direction.

